

FITNESS RANGE

Increases your strength, sculpts your body and tones your muscles.



Fit 1.0 **Fit 3.0** **Fit 5.0**

| | | | |
|----------------------------|--|--|--|
| Technology | Wired | Wired | Wireless |
| Nbr. CHANNELS | 4 | 4 | 2 (upgradeable to 4) |
| MI-Technology | NA | MI-scan | MI-scan MI-range On all channels |
| Web connection | NA | NA | Upload device history |
| Nbr. PROG. Programme offer | - Pain mgt 1 Recover 1 - Fitness 8 TOTAL: 10 | - Pain mgt 8 Recover 2 Rehab 1 Fitness 9 TOTAL: 20 | - Pain mgt 10 Recover 4 Rehab 2 Fitness 14 TOTAL: 30 |

SPORT RANGE

Builds your power, improves performance and helps you recover faster.



SP 2.0 **SP 4.0** **SP 6.0** **SP 8.0**

| | | | | |
|----------------------------|--|---|---|---|
| Technology | Wired | Wired | Wireless | Wireless |
| Nbr. CHANNELS | 4 | 4 | 4 | 4 |
| MI-Technology | MI-scan | MI-scan MI-range MI-tens | MI-scan MI-range On all channels | MI-scan MI-autorange MI-tens MI-action On all channels |
| Web connection | - | - | Upload device history | Upload device history and download objectives |
| Nbr. PROG. Programme offer | Condition 7 Pain mgt 6 Recover 3 - Fitness 4 TOTAL: 20 | Condition 10 Pain mgt 8 Recover 5 Rehab 2 Fitness 5 TOTAL: 30 | Condition 10 Pain mgt 8 Recover 5 Rehab 2 Fitness 5 TOTAL: 30 | Condition 13 Pain mgt 10 Recover 5 Rehab 2 Fitness 10 TOTAL: 40 |

SWISS TECHNOLOGY

| | |
|----------------------------|--|
| Technology | |
| Nbr. CHANNELS | |
| MI-Technology | |
| Web connection | |
| Nbr. PROG. Programme offer | |



MI, PERSONALIZED STIMULATION ADAPTED FOR EVERY PHYSIOLOGY

Equipped with unique MI (muscle intelligence) technology, Compex stimulators adapt to each of your muscles to provide the most effective electrostimulation possible, while offering greater comfort and better performance.

KREJA d.o.o.

Ljubljanska cesta 42, 1236 Trzin
Phone: +386 1 564 17 14
Mobile: +386 40 416 288
www.kreja.eu info@kreja.eu



MI, HOW IT WORKS

MI-SCAN

scans the muscle and automatically adjusts stimulator parameters to your physiology.

MI-ACTION

optimizes session effectiveness by combining voluntary contractions with electrostimulated contractions, giving you total control over your workout.

MI-TENS

makes it easier to adjust stimulation levels for pain relief programs to ensure optimal results.

MI-RANGE

indicates the optimal level of stimulation to use for Recovery and Massage programs.

MI-AUTORANGE

a brand new feature that enables your Compex module to determine the most appropriate level of stimulation automatically.

